

CB200B Charge sensitive preamplifier.

The preamplifier **CB200B** is a low noise charge sensitive preamplifier. Fast timing and small size make this preamplifier excellent module for small charged particle detectors or laboratory measurements. The preamplifier is optimized for high input capacitance (up to 1000pF). The module has bias input (up to 3 KV) and protection circuit to avoid breakdown of the input of the preamplifier circuit.



GENERAL :

Model	Charge sensitivity (Si Equivalent)	Max. Noise (keV/(Si) (Cin = 0pF)	Energy range
CB200B	60 mV/MeV	<1,2 KeV	0-130 MeV

Note: Noise characteristics see Fig.1.

PERFORMANCE

Decay time	150 us
Dynamic input capacitance:	up to 1000 pF
Noise/Input capacitance ratio:	<9 e-/pF
Integral nonlinearity:	0,1 % (without termination)
Dynamic output range :	+/- 7,5 V (without termination). +/-3 V(100 Ohms termination).
Temperature stability:	+/- 100 ppm/C.
Open loop gain:	30,000
Output resistors:	100 Ohm
Test Capacitance:	3,3 pF (+/-3%).

INPUTS/ OUTPUT

INPUT	Accepts positive or negative charge signal.
BIAS	voltage can be applied through SHV input connector. The serial resistance between input and bias connectors is 26 MegOhm.
TEST	pulse input connector is BNC type connector. Test capacitance is 3,3 pF.
POWER	input power through 3 meter screened cable from spectrometric amplifier, NIM crate power supply or portable power supply.
ENERGY	output negative or positive linear pulse. BNC connector.

POWER SUPPLY REQUIREMENTS:

The best solution is alimentation from a NIM standard power supply or special low noise linear power supplies.

Power supply pin out:

P. Voltage (V)	Current (mA)
+24	19,6
-24	10,0
+12	10,0
-12	11,6

Box dimensions: 111x80x40 mm

Pin number	
7	+24 Volt
6	-24 Volt
4	+12 Volt
9	-12 Volt
1	Ground
2	Ground

Cable length 3 m.




Fig.1. Typical noise as function of input capacitance measured with spectrometric amplifier and 2 us time constant.