

## SDIG60 single board.

The SDIG60 single module is specially designed as spectrometric amplifier. This module has gain 60 dB, zero pole, differential and integration circuits (fig.1,2). The timing of circuit has been choosed for best signal to noise resolution for specific application with CA05 preamplifier. The spectrometric chain CA05 preamplifier and SDIG60 give equivalent noise about 200 electrons. Input equivalent noise of shaper for 3 $\mu$ s (integration) and 3 $\mu$ s (differential)) signal time constant is less than 7 uV.




Fig.1 Circuit block diagram.

The module SDIG60 can be used in different applications where the signal to noise resolution is important parameter. The pole-zero cancellation with different timing to shape of input signal permit use this module board in energy spectroscopy with different detectors of particle measurement.




Fig.2. Waveform of output signal.

There are 3 trimmers on the SDIG60 board:

1. Trimmer R11 is pole-zero cancellation resistor.
2. Trimmer R7
3. adjust gain value +/- 20 %.
4. Trimmer R26 adjustment of output voltage offset.

Table 1. Pinout of connector P4.

Pin	Description
1,8	NC
2,3,5,7	GND
4	+12V
6	-12V

## PERFORMANCE

**INPUT AMPLITUDE**

**RANGE** Maximum input range 0 to  $\pm 30\text{mV}$ ; positive or negative polarity; input amplitude; protected to  $\pm 6\text{ V dc}$ ; impedance  $10\text{ k}\Omega$ , DC coupled. The range can be changed according to custom request).

**OUTPUT AMPLITUDE**

**RANGE** Unipolar non inverted. Amplitude of output signal 0 to  $\pm 10\text{ V}$  (without termination). Output impedance  $Z_o = 100\ \Omega$ . Output Offset voltage less than  $\pm 2\text{ mV}$ . The output can be terminated to  $50\text{-}\Omega$  load.

**NOISE RMS** less than  $\leq 7\ \mu\text{V}$ .

**RISE TIME**  $1,9\ \mu\text{s}$ .

**FALL TIME**  $9,8\ \mu\text{s}$ .

**NONLINEARITY**  $\leq \pm 0.05\%$  not over  $\pm 10\text{ V}$  range.  
Variation of nonlinearity in over range signal is less then  $0,1\%$ .

**TEMPERATURE**

**INSTABILITY DC** Level  $\leq \pm 25\ \mu\text{V/C}$  referred to the output (in range  $0\text{ to }50^\circ\text{C}$ ).

**COARSE GAIN** 1000.

**FINE GAIN**  $\pm 20\%$  by 25 turns trimmer.

**GAIN STABILITY** less then  $0,01\ \%/C$  in full range.

**POLE ZERO ADJ** 25-turn Trimmer. Adjustment to compensate for the preamplifier decay time constant from  $25\ \mu\text{s}$  to  $\infty$ .

**Output Offset ADJ** 14-turn Trimmer.

**TIME CONSTANT** Two stage integrate RC time constants -  $3\ \mu\text{s}$ .  
Differentiate RC time constant -  $3\ \mu\text{s}$ .

**Note:**

The RC time can be adjusted in fabric according user request.

**POWER SUPPLY REQUIREMENTS:**

Table 2.

P. Voltage (V)	Current (mA)
+12	55
-12	55

**DIMENSIONS:**

85x28x1,6 mm

**WEIGHT:**

Net 14,2 g.